

PRODUCT DESCRIPTION:

Granular Hemoglobin is obtained from desiccated bovine blood.

POTENCIAL APPLICATIONS:

It is used in preparation of bacteriological culture media as sources of nutrients for cultivation and isolation of fastidious microorganisms, especially *Neisseria* spp and *Haemophilus* spp.

PHYSICAL CHARACTERISTICS:

Granulated powder, dark reddish brown colored (almost black).

Chemical Characteristics	Specifications	Typical Value
Loss on drying	Maximum 5,00%	1,00%
pH (2% solution)	7,80 - 8,60	8,20

Growth Supporting Properties: satisfactory
 Chocolate Agar

PACKAGING	STORAGE	RETEST	CERTIFICATIONS
The dry product is packaged in polyethylene bags into reinforced fiber board drums. 25 kg 50 kg	Keep in original packaging closed, in a dry and cool place. Hygroscopic product.	4 years after its manufacturing date.	ISO 9001 SADER-SENASICA