

YEAST EXTRACT POWDER

CAT. No. 369 MI-AC-946 Version 9 Date: 14/05/2021 According to NOM-018-STPS-2015:

SECTION I. MANUFACTURER AND SUBSTANCE IDENTIFICATION

Substance name: Yeast Extract Powder.	Manufacturer name: Biotecnica Internacional S.A. de C.V.	Emergency telephone: 951 512 80 04
Other ways of identification: Saccharomyces Cerevisiae.	Address: Artículo 123 no. 104, Col. Granjas y Huertos Brenamiel, San Jacinto Amilpas, Oax. Oaxaca, México. C.P. 68285	Technical telephone: 951 512 62 55 951 512 80 04
Use and restriction: Use in culture media, fermentation processes, cell cultures and other biotechnological applications.	Issued Date: 14/08/2008	Updated: 14/05/2021

SECTION II. HAZARD IDENTIFICATION

Classification of the substance or mixture:

This substance is not classified according to the GHS. It is not classified as hazardous according to ECHA in the REACH registers (EC Regulation No. 1907/2006). No hazardous ingredients according to OSHA 29 CFR 1910.1200.

Signal word, hazards statement and precautionary statements:**Pictograms.**

It does not apply.

Signal words. (Notifications to the CLP in ECHA)

Attention: The material may cause slight irritation in prolonged or repeated contact. Do not handle it if you have known allergies or physical reaction to the yeast.

Precautionary Statement:

Prevention: P201 Obtain special instructions before use.

Prevention: P280 Wear protective gloves/protective clothing/eye protection/face protection.

Response: P301 + P312 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell.

Response: P302 + P352 IF ON SKIN: Wash with plenty of soap and water.

Response: P304 + P312 IF INHALED, call a POISON CENTER or doctor/physician if you feel unwell.

Response: P305 + P351 + P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

Storage: P402 + P404 + P235 Store in a dry and cool place. Store in a closed container.

Disposal: P501 Dispose the content or container in accordance with local/regional/national/international regulations.

Other hazards that do not contribute to the classification:

Undetermined.

SECTION III. COMPOSITION/INFORMATION ABOUT COMPONENTS

Substance name: Yeast Extract Powder.

UN number: Not applicable.

This product does not contain hazardous ingredients.

YEAST EXTRACT POWDER

CAT. No. 369 MI-AC-946 Version 9 Date: 14/05/2021 According to NOM-018-STPS-2015:

SECTION IV. FIRST AID

First aid procedures:

In case of eye contact: Remove contact lenses at the same time if you use. Rinse the eyes with plenty of water for at least 15 minutes, occasionally lifting the upper and lower eyelids. Obtain medical help if necessary.

In case of skin contact: Rinse skin with plenty of water; If irritation or adverse symptoms develop, seek medical assistance.

Ingestion: If victim is conscious and alert, rinse mouth with water as a precaution. Never give anything to an unconscious person; Get medical assistance.

Inhalation: Remove from exposure and move to fresh air immediately. If you do not breathe, give artificial respiration. If breathing is difficult, provide oxygen and seek medical help immediately.

Symptoms and health hazards:**Acute.**

Inhalation may cause respiratory irritation.

Chronic.

Prolonged contact with dust might cause irritation to respiratory tract, allergic or asthmatic reaction.

Indication of any immediate medical attention and special treatment needed:

If you have any symptoms described above, seek medical advice (show the label if possible).

SECTION V. FIRE MEASURES

Extinguishing medias:

Water, appropriate foam, carbon dioxide (CO₂), dry chemical powder.

Special hazards arising from the substance:

Any organic powder can be potentially explosive if disbursed certain concentrations. Concentrations have not been determined for these procedures.

Special fire-fighting procedures:

Data not available.

Hazardous decomposition products:

Carbon monoxide. Carbon dioxide.

SECTION VI. MEASURES TO BE TAKEN IN CASE OF SPILL OR ACCIDENTAL LEAKAGE

Personal precautions:

Property precautions must be taken to minimize the direct contact with skin, eyes and respiratory tract.

Protective equipment:

See section VIII.

Emergency procedures:

Evacuate unnecessary personal.

Environmental precautions:

Do not let the product enter sewage system.

Cleaning procedures:

Sweep the product and put it in a bag. Avoid to the dust increase, ventilate the area, and wash the place of spill after to collect all the product. Also is recommended used a vacuum system.

YEAST EXTRACT POWDER

CAT. No. 369 MI-AC-946 Version 9 Date: 14/05/2021 According to NOM-018-STPS-2015:

SECTION VII. HANDLING AND STORAGE

Cautions for safe handling:

Do not breathe dust, avoid contact with skin and eyes. Wash your hands and exposed areas with mild soap and water after using the material.

Provide good ventilation in the handling area to prevent dust buildup.

Store conditions:

Store in a clean, cool, and dry place.

Incompatibility:

Keep away from oxidizing substances, heat, flame, or sources of ignition.

SECTION VIII. EXPOSURE MEASURES/PERSONAL PROTECTION

Control parameters:

The occupational exposure limits are not available.

Technical controls:

No additional data. See Section VII.

Respiratory protection:

Follow the OSHA respirator regulations listed in 29 CFR 1910.134 or in the European Standard EN149 or equivalent. Use the NIOSH standard if the exposure limits are exceeded.

Protective gloves:

Use appropriate protective gloves to prevent exposure to the skin.

Eye protection:

Use safety glasses approved by ANSI Z87, EN166, or equivalent if a risk assessment considers it necessary.

Safety clothing:

Wear suitable protective clothing in accordance with good chemical hygiene practices.

YEAST EXTRACT POWDER

CAT. No. 369 MI-AC-946 Version 9 Date: 14/05/2021 According to NOM-018-STPS-2015:

SECTION IX. PHYSICAL AND CHEMICAL DATA

Appearance and odor: Fine yellow to light brown powder with a characteristic odor and no foreign particles.

Relative density: Not determined.

pH: 6.00 – 7.20

Solubility in water: 20 g/l

Melting Point / Freezing Point: Not determined.

Partition coefficient n-octanol / water: Not applicable.

Initial point and boiling range: Not determined.

Spontaneous ignition temperature: Not applicable.

Flammability point: Not determined.

Decomposition temperature: Not determined.

Evaporation rate: Not applicable.

Viscosity at 40 ° C (mm² / s): Not applicable.

Flammability (solid, gas): Not determined.

Average molecular weight: Data not available.

Vapor pressure (mmHg): Not applicable.

Explosive properties: Not determined.

Vapor density (air = 1): Not applicable.

Loss in drying: Maximum 5.00%.

SECTION X. STABILITY AND REACTIVITY

Reactivity:

Not established.

Chemical stability:

Stable under recommended storage conditions.

Possibility of hazardous reactions:

Not established.

Conditions to avoid:

Direct sunlight, source of ignition. Avoid humid or very warm places.

Incompatible materials:

Strong acids. Strong bases.

Hazardous decomposition products:

Fumes. Carbon monoxide. Carbon dioxide.

YEAST EXTRACT POWDER

CAT. No. 369 MI-AC-946 Version 9 Date: 14/05/2021 According to NOM-018-STPS-2015:

SECTION XI. TOXICOLOGICAL INFORMATION

Signs and symptoms of exposure:

It may cause irritation of the respiratory tract, skin and eyes. Risks are unlikely with proper handling, respecting safety practices.

Acute toxicity estimates:

Not Classified.

Specific toxicity in certain organs:

Not Classified.

Pathways of entry (See section II):

Skin: Yes

Eyes: Yes

Inhalation: Yes

Ingestion: Yes

Carcinogenicity: Not listed in the NTP report

IARC monographs: Not listed.

OSHA 1910.1002: Not listed.

Interactive effects: Data not available.

SECTION XII. ECOTOXICOLOGICAL INFORMATION

Toxicity: This product is not hazardous to the environment according to US regulations.

Persistence and Biodegradability: Information not available.

Bioaccumulation potential: Information not available.

Mobility on the ground: Information not available.

PBT and vPvB: The substance does not meet the criteria to be identified as PBT substance or vPvB according to Annex XIII of the REACH regulation.

Other adverse effects: Has not been evaluated.

SECTION XIII. INFORMATION RELATING TO THE DISPOSAL OF THE PRODUCTS

Contact a licensed professional waste disposal service to dispose of this material. Avoid release to the environment or enter the sewage system. Check all local environmental regulations, estate and federal. See section VIII.

YEAST EXTRACT POWDER

CAT. No. 369 MI-AC-946 Version 9 Date: 14/05/2021 According to NOM-018-STPS-2015:

SECTION XIV. INFORMATION CONCERNING TRANSPORTATION

UN Number:

Not Applicable.

Official Transport Designation of UN:

Not Applicable.

Hazards in transportation:

Not Applicable.

IATA:

Non-Hazardous for Air Transportation.

ADR/RID:

Non-Hazardous for Road Transportation.

IMDG:

Non-Hazardous for Sea Transportation.

Environmental Risks:

Have not been thoroughly investigated, but currently none have been identified

Special Precautions for Users:

Not Classified.

Transport in Bulk According to Annex II of Marpol 73/78 and the IBC Code:

Not Classified.

SECTION XV. REGULATORY INFORMATION

Not applicable.

SECTION XVI. OTHER INFORMATION

The information is considered correct but is not exhaustive and will be used only as guidance, which is based on the current knowledge of the chemical or mixture and is applicable to the appropriate safety precautions for the product.

YEAST EXTRACT POWDER

CAT. No. 369 MI-AC-946 Version 9 Date: 14/05/2021 According to NOM-018-STPS-2015:

Abbreviations or acronyms:

EC: European Inventory of Existing Chemical Substances.

UN Number: Number of four digits used to identify hazardous substances or materials.

NIOSH: National Institute of Occupational Safety and Health of the United States.

FFP1: Level of protection of a mask according to the European standard EN 14.

IN: European Standard.

ANSI: American National Standards Institute.

NTP: National Toxicology Program of the United States.

IARC: International Agency for Research on Cancer.

OSHA: Occupational Safety and Health Administration in the United States.

PBT: Persistent, Bioaccumulative and Toxic.

vPvB: Very persistent and very Bioaccumulable.

REACH: Registration, Evaluation, and Authorization and Restriction of Chemical substances of the European community.

IATA: International Air Transport Association.

ADR / RID: European Agreement on the International Carriage of Dangerous Goods by Road / European Agreement on the International Carriage of Dangerous Goods by Rail.

IMDG: International Maritime Dangerous Goods Code.

MARPOL: International Convention for the Prevention of Pollution from Ships.

CIQ: International Code of Chemotherapy.

ECHA: European Chemicals Agency for Chemical Substances and Mixtures.

TSCA: Toxic Substances Control Act Chemical Substance Inventory contains all existing chemical substances manufactured, processed, or imported into the United States.

GHS: Globally Harmonized System.